



**We make your Lab
a safer Place!**

World Premiere:

Safety Waste Cap

for small Neck Threads



Safety Waste Cap for small Neck Threads

Discover our most flexible Waste System!

✓ Fits smallest thread sizes

✓ Use empty solvent bottles for HPLC waste

Turn your empty solvent supply bottles into a safe collection container for HPLC waste!

✓ Suitable for MERCK and Wheaton bottles

Including thread sizes GL/S 40 for MERCK and 38-430 for Wheaton bottles.

✓ Easily customizable

- Connects every HPLC system safely to the waste container.
- Blind plugs included for every connection - for a safely closed system.

✓ Workspace Safety

- Connection for exhaust filter - keeps your lab air clean and protects your health.

3 Tubing connectors

- For 5.0 - 11.5 ID tubing
- Including blind plugs for connectors not in use

4 Capillary connectors

- Standard connections for HPLC capillaries with 1.6mm, 2.3mm and 3.2mm OD
- Including all 3 fitting sizes for each connection



Easily expandable

- Put additional HPLC manifolds on top do double the number of connections!

Threads

• GL 28	• GL 38
• GL 32	• 38-430
	• S 40


Safety Waste Cap for small Neck Thread Bottles

Ordering Information

A 350 028



B 350 032




C 350 038




D 350 039



E 350 040



F 160 532



Pic.	Part no.	Description	Thread size	Connections for 1.6/2.3/3.2mm OD tubing	Connections for 5.0-11.5 ID tubing	Unit
A	350 028	Safety Waste Cap Mini-LISA, GL 28	GL 45	4x	3x	1
B	350 032	Safety Waste Cap Mini-LISA, GL 32	GL 45	4x	3x	1
C	350 038	Safety Waste Cap Mini-LISA, 38-430	GL 45	4x	3x	1
D	350 039	Safety Waste Cap Mini-LISA, GL 38	GL 45	4x	3x	1
E	350 040	Safety Waste Cap Mini-LISA, GL/S 40	GL 45	4x	3x	1
F	160 532	HPLC Manifold PTFE, GL 14	GL 14	4x	3x	

Safety Waste Cap for small Neck Thread Bottles

Consumables & Accessories

A 107 063**B** 107 059**C** 107 061**D** 160 502**E** 117 808**F** 160 506**G** 107 620**H** 490 335**I** 490 336**J** 407 986

Pic.	Part no.	Description	Diameter	Thread size	Material	Color	Unit
A	107 063	PFA fitting with integrated ferrule	3.2 mm ID	UNF 1/4" 28G	PFA	Blue	5
B	107 059	PFA fitting with integrated ferrule	2.3 mm ID	UNF 1/4" 28G	PFA	Violet	5
C	107 061	PFA fitting with integrated ferrule	1.6 mm ID	UNF 1/4" 28G	PFA	Green	5
D	160 502	Blind plug, for capillary connector	-	UNF 1/4" 28G	PFA	Natural	5
E	117 808	Tube connector, curved	5.0-11.5 mm OD	NPT 1/8"	PP	Blue	1
F	160 506	Blind plug for tube connector	-	NPT 1/8"	PP	Blue	1
G	107 620	Blind plug, for exhaust filter connector	-	GL 14	PTFE	White	1
H	490 335	Exhaust filter S, service life 3 months each	-	GL 14	PP	Red	4
I	490 336	Exhaust filter M, service life 6 months each	-	GL 14	PP	Red	2
J	407 986	Exhaust filter L, service life 12 months each	-	GL 14	PP	Red	1



SCAT - Europe GmbH


Waldecker Straße 7
64546 Mörfelden-Walldorf
Deutschland

 **+49 (0) 61 05 / 30 55 86 - 0**

 **+49 (0) 61 05 / 30 55 86 - 99**

 **info@scat-europe.com**

 **www.scat-europe.com**



SCAT Lab Safety

1704 East Blvd., Suite 101
Charlotte, NC 28203
USA

 **+1-800-454-1326**

 **info@scatlabsafety.com**

 **www.scatlabsafety.com**

